LIVSANE

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EC REP

MedPath GmbH
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

Shenzhen Jumper Medical Equipment Co., Ltd.
D Building, No. 71, Xintian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China 518103
Tel:+86-755-26696279
E-mail: info@jumper-medical.com
Website: www.jumpermmed.com
www.jumper-medical.com

CE 0598

SI

EU: PXG Pharma GmbH,
Pfingstweidstrasse 10-12,
68199 Mannheim,
Germany.
www.pxgpharma.com

LIVSANE

EAR AND FOREHEAD THERMOMETER

Instruction Manual

Modell: JPD-FR301

EN

Contents

Unpacking Check 57
Package Contents 57
Safety Precautions 58
Warning 59
Symbols 60
Body Temperature Basics 60
Product Description 61
Features 62
Product Structure 63
Display Description 63
Sounds and Backlight Color Instructions 64
Display and Operating Instructions 65
Measuring Ear Temperature 70
Measuring Forehead Temperature 70
Measuring Object Temperature 71
Replacing Batteries 72
Cleaning and Disinfection 73
Maintenance 74
Troubleshooting 74
Specifications 75
Security Class 76
Storage and Transportation 76
EMC Information-Guidance and Manufacture's Declaration 76
Warranty and After-Sale Service 79

Product Information

Product Name: Ear and Forehead Thermometer

Model: JPD-FR301

Manufacturer: Shenzhen Jumper Medical Equipment Co., Ltd

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China. 518103

Copyright

Copyright© Jumper Medical. All rights reserved.

Statement

Shenzhen Jumper Medical Equipment Co., Ltd. owns and reserves all of the rights comprised in the copyright of this document. No part of this document may be changed, excerpted, copied, reproduced, or imitated in any form or by any means without the prior consent of Shenzhen Jumper Medical Equipment Co., Ltd.

All statements, information, and recommendations in this document are provided "AS IS" without warranties, guarantees or representations of any kind, either express or implied. The information in this document is subject to change without notice. Shenzhen Jumper Medical Equipment Co., Ltd. reserves the right of final interpretation of this document.

Introduction

Thank you for purchasing this Ear and Forehead Thermometer. Please read the User Manual carefully to make sure safe and proper use of this thermometer.

Please read and fully understand the Safety Precautions before use.

Keep the Instruction Manual with this thermometer for future reference.

Unpacking Check

Please open the package carefully before use, check whether all accessories are available or not and whether any component is damaged during transportation, and perform installation and operation following this user manual. In case of any damage or operation problem, please contact the dealer or contact Jumper directly. You will need the following information when making your claim: device model, serial number, purchase date, and your contact information and address.

Package Contents

N°.	Name	Quantity
1	Ear and Forehead Thermometer	1
2	Pouch	1
3	Battery (AAA)	2
4	Instruction Manual	1

Safety Precautions

Read the following precautions carefully before using the hermometer.

⚠ Attention
• Take care of the temperature probe lens, which is fragile.
• No service /maintenance while the equipment is in use.
• Dispose used batteries with care. To protect the environment, you are recommended to send the used batteries to a designated collection point.
• The thermometer doesn't need recalibration.
• - Make sure that the thermometer is not exposed to lint, dust, light (including sunlight), etc. - Please note the effects of degraded sensors that can degrade performance or cause other problems. - Make sure that the thermometer is not exposed to pets or pests.
• If the thermometer is soiled or its infrared optical components is damaged, please stop using it.
• The lay operator or lay responsible organization should contact the manufacturer or manufacturer's representative on the following issues: - assistance in setting up, using, or maintaining the equipment or system when needed, or to report unexpected operation or events.
• The IR thermometer is identified as intended for consumer use. The patient is an intended operator.
• The effect that the following actions could be unsafe as applicable: use of accessories, detachable parts, and materials not described in the instructions for use or modification of the equipment.
• The hazards that can result from unauthorized modification of the equipment.
• The manufacture can provide the circuit diagram, component part list, description and calibration instructions to assist service personal for parts repair.
• Do not subject the thermometer to vibration or impact.
• Do not take body temperature readings within 20 minutes after you do physical exercises or get excited.
• Do not use the thermometer on newborns or for continuous temperature monitoring purposes.
• Do not use the thermometer for purposes that are not specified in this User Manual. Follow the instructions in the "Measurement Process" chapter and carefully operate the thermometer when measuring children's temperature.
• Do not immerse the thermometer in water or other liquid. Clean and disinfect the thermometer as described in the "Cleaning and Disinfection" chapter.
• Do not touch the tip of the temperature probe, on which a precise temperature sensor resides.

• Keep the temperature probe clean to make sure accurate readings.
• Before measuring the temperature from the ear canal, remove any earwax.
• The ambient temperature must not be extremely high or low. To make sure accurate readings, keep the thermometer under room temperature for more than 30 minutes before use.
• Do not use the thermometer under an ambient temperature higher than 40°C (104°F) or lower than 10°C (50°F), which is beyond the operating temperature range of the thermometer.
• **Risk of pollution!** The user is recommended to send the overdue thermometer to local garbage disposal site or send it back to us.
• 2 AAA batteries of 1.5V are the only replaceable accessories of the thermometer. Please do not use the batteries of other voltages or specifications.

Warning

⚠ **Warning**
• Do not force the temperature probe of the thermometer into an ear canal. Otherwise, the ear canal may get injured.
• Do not use the thermometer if the ear is infected with otitis or suppuration.
• Keep the thermometer out of the reach of children.
• The result may be inaccurate if you use the overdue thermometer.
• The thermometer is not intended to diagnose or treat any health problem or disease. The measurement results are for reference only.
• It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results. For such purposes, please consult a doctor.
⃠ Do not charge an alkaline dry-cell battery or throw it in fire. Otherwise, the battery may explode.
⃠ Do not disassemble the thermometer or attempt to repair it. Otherwise, the thermometer may be damaged permanently.
⃠ Do not take temperature measurements on body parts other than forehead and ears. Otherwise, the temperature readings may be inaccurate.
⃠ During measurement, do not use a mobile phone or any other device that may cause electromagnetic interference.
⃠ Do not use the thermometer in an environment where flammable anesthetic mixture with air or with oxygen or nitrous oxide is available.
⃠ Do not open the battery cover when using the thermometer.

Symbols

Symbol	Description
	Type BF applied part.
	Attention must be paid.
	The action is prohibited.
	Information about the manufacturer.
	Date of manufacture.
	Consult the instructions for use.
MD	Medical Device
CE 0598	This product complies with the Regulation (EU) 2017/745 requirements.
	Waste electrical materials should be sent to a dedicated collection point for recycling.Must contact its local authorities to determine the proper method of disposal of potentially bio hazardous parts and ACCESSORIES.
IP22	Degree of protection against the Ingress of water.
	Humidity
	Temperature limitation
	Atmospheric Pressure
 Avertissement	A personal injury or damage to the thermometer may occur if the thermometer is not used correctly.
 Attention	Inaccurate reading or damage to the thermometer may occur if the thermometer is not used correctly.
UDI	Unique Device Identifier
EC REP	Authorized European Representative
SN	Serial number
LOT	Batch code

Body Temperature Basics

• The normal body temperature is a range.
• The normal range varies from person to person and can fluctuate throughout the day.
• The normal range also varies by body site. Therefore, measurements from different sites should not be compared directly.


EN

To determine if an individual is experiencing an elevated body temperature and/or having a fever, it is critical to know the individual's normal body temperature when he/she is well. Take multiple readings to obtain the normal body temperature range and note the specific body site measured, for example: forehead or eardrum temperature.

Body Site	Normal Temperature Range
Forehead	34,7°C-37,3°C (94,5°F- 99,1°F)
Eardrum	35,8°C-38,0°C (96,4°F-100,4°F)
Mouth	35,5°C-37,5°C (95,9°F- 99,5°F)
Armpit	34,7°C-37,3°C (94,5°F- 99,1°F)
Rectal	36,6°C-38,0°C (97,9°F-100,4°F)

The normal body temperature range varies slightly with age and gender. Generally, newborns or children have higher body temperatures than adults, and adults have higher body temperatures than the elderly. Women's body temperatures are approximately 0.3°C (0.5 °F) higher than men's.

Variation in body temperature

Normal body temperature fluctuates throughout the day and is also affected by external factors.The body temperature of an individual is the lowest between 2:00 a.m. and 4:00 a.m. and the highest between 2:00 p.m. and 8:00 p.m. An individual's body temperature typically changes by less than 1°C (1.8°F) each day.

Product Description

1. Overview


Ear and Forehead Thermometer JPD-FR301 measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can quickly get measurement results after positioning properly the temperature probe in the ear canal or forehead.It can also be used to measure the surface temperature of object.(such as milk and water).

2. Structure

The thermometer consists of a shell, an LED, a measure button, a beeper, an infrared temperature sensor, and a Microprocessor.

3. Operating principle

The Ear and Forehead Thermometer sensor collects infrared energy emitted by the eardrum or the skin surface. After being focused by a lens, the energy is converted into a temperature reading by the thermopiles and measurement circuits.

4. Intended use

The Ear and Forehead Thermometer takes human body temperature via the eardrum or forehead for both professional and home use.

5. Contraindications

The device has no side-effects if administered correctly and residual risk is acceptable.

Features

01. Good safety

- Passive infrared receiving technology.

02. Easy operation

- Ergonomic design
- One button measurement

03. Quick measurement

- 1-second measurement

04. High accuracy

- Advanced infrared temperature sensor, with high sensitivity
- High accuracy with automatic temperature calibration

05. Diverse functions

- R20 readings recall
- Fever alert
- Switching between °C and °F
- Automatic power-off, power saving


06. Intended patient population

- Forehead mode applies to all age groups
- Ear mode applies to people above 3 months

07. Child mode

- This mode is recommended for people under 12 years of age.

Product Structure


1. LED display screen
2. Power button / Measure button
3. Memory button / Sound switch
4. Mode button (Adult mode / Child mode/ Object)
5. Probe (take off the cover when measuring the ear temperature)(applied part)
6. Unit switch button (˚C / ˚F)
7. Probe cover (Automatic switch among object, forehead and ear modes, by removing the probe cover to measure ear temperature and placing back the probe cover for object and forehead temperature measurement.)
8. Battery cover

EN

Display Description

1. Object temperature mode
2. Forehead temperature mode
3. Child mode
4. Ear temperature mode
5. Mute / un-mute
6. Memory recall
7. Temperature unit (˚C)
8. Temperature unit (˚F)
9. Low battery
10 Temperature value

Sounds Instructions

Range	Sounds
Forehead temperature (Adult / Child)	
22.0°C-37.5°C/71.6°F-99.5°F	A long beep
37.6°C-43°C/99.6°F-109.4°F	3 short double beeps
Ear temperature (Adult / Child)	
34.0°C-37.5°C/93.2°F-99.5°F	A long beep
37.6°C-43.0°C/99.6°F-109.4°F	3 short double beeps
Object temperature	
0°C-100°C/32.0°F-212°F	A long beep

Note: When the forehead temperature is between 22.0°C/71.6°F and 37.5°C/99.5°F, the ear temperature is between 34.0°C/93.2°F and 37.5°C/99.5°F, there will be a long beep. This indicates that your body temperature is normal.


When the forehead and ear temperature is between 37.6°C/99.6°F and 43.0°C/109.4°F, there will be 3 short double beeps and the temperature value will flash on the screen. This indicates that your body temperature is a high. You may have a fever. Please consult your doctor if you are not sure.



EN

Display and Operating Instructions


Screen Display	Operating Instructions Displayed State	Sound
Measuring Ear temperature (Adult / Child)		
  Ear temperature for adult  Ear temperature for child	Take off the probe cover, press and release the **Power button** for 1 second to power on the thermometer. The symbol " " is displayed on the screen. Switch to the adult or child measurement mode by pressing the **Mode button** according to your measurement needs. Insert the temperature probe into a proper position in the ear canal. Press the **Measure button** to start a measurement.	See the table in the "Sounds instructions" section.



EN

Screen Display	Operating Instructions Displayed State	Sound
Measuring Forehead temperature (Adult / Child)		
 Forehead temperature for adulte  Forehead temperature for child	Put the cover on the probe, press and release the **Power button** for 1 second to power on the thermometer The "Head" symbol is displayed on the screen. Switch to the adult or child measurement mode by pressing the **Mode button** according to your measurement needs. Point the thermometer to the center of the forehead, about"1-3cm " away from the skin surface. Press and release the **Measure button**. The temperature reading will be displayed on the screen.	See the table in the "Sounds instructions" section.
Measuring Object temperature		
	Put the cover on the probe, switch to the object measurement mode by pressing the **Mode button**. The thermometer enters the **Object mode**. The "House" symbol is displayed on the screen. Point the thermometer to the center of the object. Press and release the **Measure button**. The temperature reading will be displayed on the screen.	See the table in the "Sounds instructions" section.
Out of the measuring range display		
 (for reference only)	In Ear mode, a temperature reading of more than 43.0°C (109.4°F) In Forehead mode, a temperature reading of more than 43.0°C (109.4°F) In Object mode, a temperature reading of more than 100°C (212.0°F)	A long beep


Screen Display	Operating Instructions Displayed State	Sound
 (for reference only)	In Ear mode, a temperature reading of less than 34.0°C (93.2.0°F) In Forehead mode, a temperature reading if less than 22.0°C (71.6°F) In Object mode, a temperature reading of less than 0°C (32.0°F)	A long beep
Recall 20 memories		
	In a power-on state, press the **Memory butto**n to enter the memory mode. When the **Memory button** is released, **01** will be shown, followed by the recorded reading. Press the **Memory button** again for the next recorded data, **02** will be shown, followed by the recorded reading. A maximum of 20 temperature readings can be recalled. When the maximum number of records is exceeded, the earliest memory data will be overwritten. **Note:** **01** means the latest data.	Silent
No memory data / Clear memory data		
	The display is as shown , when there is no more data checked while recalling memories. Remove 2 dry batteries and after 10 seconds re-install the power to clear all memory data.	When the power is turned back on, there are a long beep

EN

Screen Display	Operating Instructions Displayed State	Sound
Switching between mute and un-mute		
	In the power-on state, press and hold the sound switch button for about 2 seconds to switch the sound on or off. When the sound is turned on, it will beep once and the mute symbol will be displayed when the sound is off. The 🔇 symbol shows up in Mute mode and disappears in Un-mute mode.	When the sound is turned on, it sounds a long beep
Switching between object temperature and body temperature		
	When putting on the probe cover, press the **Mode button** to switch between object and forehead modes. Ear mode (adult/child) will be activated by removing the probe cover.	Silent


Screen Display	Operating Instructions Displayed State	Sound
Switching between °F/°C		
	Remove the battery cover, and toggle the unit switch button in the battery compartment to switch the °C/°F unit.	Silent
Error information & low battery		
	The ambient temperature is higher than 40.0°C (104.0°F) or lower than 10.0°C (50.0°F).	A long beep
	An error occurs when data is being read from or written to the memory, or the temperature correction is not complete.	A long beep
	When the battery voltage is lower than 2.5V ± 0.1V, the low battery symbol will appear on the display. Please replace the batteries.	Silent

Measuring Ear Temperature

1. When using the thermometer for the first time, install the battery according to the stated polarities.
2. Take the probe cover off from the thermometer before measuring the ear temperature.
3. Press the **Power button** to power on the thermometer.
4. The "Ear" symbol is displayed on the screen.
5. Switch between "adult" and "child" modes by pressing Mode button symbol shows up in "child" mode.
6. Insert the temperature probe into the ear canal.
7. Press and release the **Measure button**. The ear temperature reading will be displayed on the screen instantly.

Note: Children under 1 year: Pull the ear straight back.
Children aged 1 year to adult: Pull the ear up and back.

⚠ Do not force the thermometer into the ear canal. Otherwise, the ear canal may get injured.


⚠ When taking the temperature on an adult, gently pull the ear up and back to make sure the ear canal is straight, so that the temperature probe can receive an infrared ray from the eardrum.


⚠ Be careful when taking temperature on a child, whose ear canal is small.

Measuring Forehead Temperature

1. Put the cover on the probe of thermometer.
2. Press the **Power button** to power on the thermometer.
3. Press the **Mode button**, the thermometer enters the **Forehead mode**. The " " symbol is displayed on the screen.
4. Switch between "adult" and "child" modes by pressing **Mode button**. symbol shows up in "child" mode.
5. Point the thermometer probe to the center of the forehead, about " 1-3cm "away from the skin surface.

EN


6. Press and release the **Measure button** for 1 second. The temperature reading will be displayed on the screen instantly.
7. If no activity is detected, the thermometer will power off automatically in 10 seconds.

Measuring Object Temperature

1. Put the cover on the probe of thermometer.
2. Press the **Power button** to power on the thermometer.
3. Press the **Mode button**, the thermometer enters the **Object mode**. The " " symbol is displayed on the screen.
4. Point the thermometer probe to the center of the object, about "1-3cm(0,39 - 1,18 inches)" away from the object surface.
5. Press and release the **Measure button** for 1 second. The temperature reading will be displayed on the screen instantly.
6. If no activity is detected, the thermometer will power off automatically in 10 seconds.

After a measurement

(1) After each measurement, you can enter the recall mode and query earlier temperature readings. For more details, see "Recall 20 memories" in the preceding table.

(2) After each measurement, clean the temperature probe with a soft cloth, and put the thermometer in a dry and well-ventilated place.

Notes:

(1) The thermometer is suitable for an indoor environment without strong air convection (for example, winds from a fan, an air-conditioner, or a heater) between the thermometer and the person.

(2) Make sure that the ear canal is clean and dry before starting a measurement. It is recommended to clean the ear canal with a cotton swab if any dirt exists. Otherwise, the temperature probe may be polluted and temperature readings may be inaccurate.

(3) Do not hold the thermometer for a long time, because it is sensitive to the ambient temperature.

(4) Make sure the sensor is free of foreign matters before use.


(5) Make sure the forehead has no sweat and no hair covered before measuring the forehead temperature; otherwise, the result could be incorrect.

(6) No intense emotion or strenuous exercises before measuring.

(7) After measuring the data once, you should wait for the unit symbol (˚C / ˚F) flashing and measure the next data.

Replacing Batteries

1. Slide the battery cover off along the marked direction and take it off.
2. Insert the two AAA batteries into the compartment according to the stated polarities.

⚠ Make sure that the batteries are installed correctly. Otherwise, the thermometer may be damaged.

⚠ If the low-battery symbol is displayed on the screen, replace the batteries.

⚠ Batteries of same type should be used. Dispose the used batteries in accordance with the local environmental policies.

⚠ The thermometer is shipped with batteries. First open the battery cover, then remove the insulating piece.

EN


Cleaning and Disinfection

Cleaning

Recommended detergents:

* Medical detergents
* Home use mild detergents

Cleaning steps:

(1) Take the batteries out before cleaning.
(2) Clean the temperature probe with a soft cloth. Clean the lens of the temperature probe with a cotton swab.
(3) Wipe the thermometer body with a slightly damp soft cloth.

⚠ Keep water out off the lens during the cleaning process. Otherwise, the lens may be damaged.

⚠ The lens may be scratched if it is cleaned with a hard object, which might result in inaccurate readings.

⚠ Do not clean the thermometer with corrosive cleansers. During the cleaning process, do not immerse any part of the thermometer into liquid, or allow liquid to penetrate the thermometer.

Disinfection

Recommended disinfectants:

*Isopropyl alcohol solution (concentration: 70%)
*Medicinal alcohol (concentration: 75%)
*Sodium hypochlorite solution (concentration: 3%)

Disinfecting steps: :

1. Wet the clean soft cloth with a small quantity of disinfectant, wipe the thermometer and quickly dry it.
2. Disinfect the thermometer body and the area around the temperature probe with a cloth slightly moistened with 75% medical alcohol.

⚠ Do not use hot steam or ultraviolet radiation for disinfection. Otherwise, the thermometer may be damaged or quickly aged.

⚠ It is recommended to disinfect the thermometer before and after each use. The disinfection time is completed within 1 minute, and the number of repetitions per disinfection is not more than 2 times.

⚠ Clean and disinfect the thermometer under the temperature of +10°C~+40°C (50°F-104°F), the relative humidity of 15%~85%RH (no condensation) and the barometric pressure of 86kPa~106kPa.

EN

Maintenance

Preventive inspection & maintenance period

1. Ensure the safety of thermometer, and check whether it has potential safety hazards in normal use each week, e.g. whether the lens is broken, the shell has cracks and the sensing head is polluted. Do not use the thermometer with potential safety hazard. Clean the thermometer if not used for a long time.
2. After each use, clean the temperature probe as described in the "Cleaning and Disinfection" chapter.
3. Store the thermometer in a dry, dust-free, and well-ventilated place. Make sure that the thermometer is not exposed to sunlight. Make sure that the storage and transportation environments meet the requirements.
4. Check regularly whether safety risks exist.
5. Remove the batteries if the thermometer will not be used for more than two months.

Troubleshooting

Problem	Possible Cause	Solution
The thermometer fails to power on	Low battery	Change the batteries.
	Polarities of the batteries are reversed.	Make sure that the batteries are installed correctly.
	The thermometer is damaged.	Contact the manufacturer.
"Er1" is displayed	The ambient temperature is lower than 10°C (50.0°F) or higher than 40°C (104°F).	Take a measurement under an ambient temperature between 10°C (50.0°F) and 40°C (104°F).
The temperature reading is lower than the typical body temperature range.	The lens of the temperature probe is dirty.	Clean the lens using a cotton swab.
	The thermometer probe is not aligned to the eardrum.	Reposition the thermometer probe so that it is aligned to the eardrum.
	The thermometer is used within 30 minutes after being taken from a cold environment.	Wait for more than 30 minutes after the thermometer is moved into the measurement environment.
The temperature reading is higher than the typical body temperature range.	The temperature probe is damaged.	Contact the manufacturer.

Specifications

Product Name	Ear and Forehead Thermometer
Product Model	JPD-FR301
Power Supply Mode	Internal power supply
Operating Voltage	DC 3V
Battery Model	AAA x 2
Operating Mode	Continuous operating
Display	LED screen
Measure time	About 1 second
Latency Time	About 3 seconds
Measuring Range	Forehead mode:22.0℃-43.0℃ (71.6℉-109.4℉) Ear mode:34.0℃-43.0℃ (93.2℉-109.4℉) Object mode:0.0℃-100.0℃ (32.0℉-212.0℉)
Accuracy (Laboratory)	Forehead mode: ±0.2℃ (36.0℃-39.0℃) ; ±0.3℃ (22.0℃-36.0℃ / 39.0℃ ~43.0℃); Ear mode: ±0.2℃ (36.0℃-39.0℃) ; ±0.3℃ (34.0℃-36.0℃ / 39.0℃~43.0℃); Object mode: ±1.0℃/±2.0℉
Accuracy (Clinical)	±0,3℃ (±0,6℉)
Measuring Site	Ear canal, Forehead (keep distance 1~3 cm from forehead)
Reference Body Site	Armpit
Resolution	0,1℃ (0,1℉)
Mode of Operation	Adjusted mode
Memory	20 temperature readings
Low-battery Alert	The low-battery symbol is displayed if the power voltage is lower than 2.5 V±0.1V
Automatic Power-off	The thermometer automatically powers off if it is not used in 10±1 seconds.
Outer dimensions (mm)	145.5*43.4*45.4mm
Weight (g)	Thermometer (with batteries): 94g
Manufacturing Date	See the label
Software version	1.0.0
Service Life	2 years
Battery Life	Alkaline dry battery for around 20,000 measurements
Operating Environment	Temperature: 10℃~ 40℃ (50℉-104℉)
	Humidity: 15%-95% RH, non-condensing
	Atmospheric pressure: 86-106 kPa

EN

Storage and Transportation	Temperature: -20°C to 50°C (-4°F-122°F)
	Humidity:15- 95% RH, non-condensing
	Atmospheric pressure: 50-106 kPa

The Ear and Forehead Thermometer has been tested and conforms to the standard ASTM E1965-98. ASTM laboratory accuracy requirements in the display range of 96.8°F to102.2°F (36°C-39°C) for ear canal IR thermometers is ±0.4°F (±0.2°C). Note that for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is ±0.2°F (±0.1°C).

Security Class

Type of protection against electric shock: internally powered equipment.

Degree of protection against electric shock: Type BF applied part.

- Degree of protection against ingress of water: IP22
- Safety degree of using in flammable anesthetic gas blending with air, oxygen or nitrous oxide: Non-AP/APG
- No application parts of the thermometer prevents defibrillation charge effect.
- No application parts of the thermometer output signal.
- The thermometer is impermanent installed device.

Storage and Transportation

The thermometer can be transported using general transportation tools. Severe vibration, shock, or rain must be avoided during transportation.

The thermometer must be packaged and then stored in a well-ventilated room without corrosive gas.

EMC Information-Guidance and Manufacture's Declaration

1* WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally."

2* WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

3* WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result."

Table 1

declaration - electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Table 2

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°C, 45°C, 90°C, 135°C, 180°C, 225°C, 270°C and 315°C 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°C 0 % UT; 250/300 cycles	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

EN

Table 3

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10 V/m

Table 4

declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maximum	Niveau d'immunité	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MH	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m
Note* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.					

Warranty and After-Sale Service

The device is under warranty for two years from the date of purchase.

The batteries, the packaging, and any damage caused by improper use are not covered by the warranty.

Excluding the following user-caused failures:

1. Failure resulting from unauthorized disassembly and modification.
2. Failure resulting from an unexpected dropping during application or transportation.
3. Failure resulting from not following the instructions in the user's manual.

After-sale service unit:
Shenzhen Jumper Medical Equipment Co., Ltd.
Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen,Guangdong,China 518103
Tel: +86-755-26696279
Fax: +86-755-26852025
E-mail: info@jumper-medical.com
Website: www.jumpermmed.com
www.jumper-medical.com
Postal Code: 518103

EC REP

MedPath GmbH
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

Shenzhen Jumper Medical Equipment Co., Ltd.
D Building, No. 71, Xintian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China 518103
Tel:+86-755-26696279
E-mail: info@jumper-medical.com
Website: www.jumpermmed.com
www.jumper-medical.com

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Copyright© Jumper Medical. Tutti i diritti sono riservati.
Copyright© Jumper Medical. Vse pravice pridržane.